

- ✓ - Indicates «Combination approved by Medacta International»
- ✗ - Indicates «Combination NOT approved by Medacta International»

*This represents functional compatibility and does not necessary indicate regulatory approval status for marketing of product combinations.
For further product regulatory information please contact your Medacta representative.*

		Stems													
		AMIS-K (12/14 taper)	AMISstem C (12/14 taper)	AMISstem H (12/14 taper)	AMISstem H Collared (12/14 taper)	AMISstem H Proximal Coating (12/14 taper)	MasterLoc (12/14 taper)	MiniMAX (12/14 taper)	M-Vizion Femoral Revision System (12/14 taper)	Native (10/12 taper)	Quadra-C (12/14 taper)	Quadra-H (12/14 taper)	Quadra-S (12/14 taper)	Quadra-R (12/14 taper)	X-Acta (12/14 taper)
Heads	Stainless Steel (10/12 taper)	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗
	Stainless Steel (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	CoCr (10/12 taper)	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✗	✗	✗	✗
	CoCr (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Mectacer BIOLOX <i>delta</i> (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Mectacer BIOLOX OPTION (ball head + sleeve 12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Ceramtec BIOLOX <i>delta</i> (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Ceramtec BIOLOX OPTION (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Endo Head (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓

- ✓ - Indicates «Combination approved by Medacta International»
- ✗ - Indicates «Combination NOT approved by Medacta International»

		Liners							UHMWPE Cups		
		Versafitcup Double Mobility UHMWPE / DM Liner UHMWPE	Versafitcup Double Mobility HC UHMWPE / DM Liner HC UHMWPE	Versafitcup CC Flat and Hooded UHMWPE	Versafitcup CC Flat and Hooded HC UHMWPE	Mectacer BIOLOX <i>delta</i>	Ceramtec BIOLOX <i>delta</i>	Impact Flat and Hooded HC UHMWPE	Apricot Cemented HC UHMWPE flat cup	Native Cemented HC UHMWPE Cup	Medacta Bipolar Head
Heads	Stainless Steel (10/12 taper)	✓	✓	✓	✓	✗	✗	✓	✓	✓	✗
	Stainless Steel (12/14 taper)	✓	✓	✓	✓	✗	✗	✓	✓	✓	✗
	CoCr (10/12 taper)	✓	✓	✓	✓	✗	✗	✓	✓	✓	✗
	CoCr (12/14 taper)	✓	✓	✓	✓	✗	✗	✓	✓	✓	✗
	Mectacer BIOLOX <i>delta</i> (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗
	Mectacer BIOLOX OPTION (ball head + sleeve 12/14 taper)	iii	iii	✓	✓	✓	✓	✓	✓	✗	✗
	Ceramtec BIOLOX <i>delta</i> (12/14 taper)	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗
	Ceramtec BIOLOX OPTION (12/14 taper)	iii	iii	✓	✓	✓	✓	✓	✓	✗	✗
	Endo Head (12/14 taper)	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗

iii - The internal sleeve of BIOLOX OPTION of ø28 mm, size XL, and sometimes size L, do not completely cover the femoral stem threaded taper. This may cause a premature wear of the double mobility liner. We recommend to follow the combinations indicated in the table below "Double Mobility Liners - BIOLOX OPTION Combinations".

DOUBLE MOBILITY LINERS - BIOLOX OPTION COMBINATIONS

		Stems												
		AMIS-K (12/14 taper)	AMISstem C (12/14 taper)	AMISstem H (12/14 taper)	AMISstem H Collared (12/14 taper)	AMISstem H Proximal Coating (12/14 taper)	MasterLoc (12/14 taper)	MiniMAX (12/14 taper)	M-Vizion (12/14 taper)	Quadra-C (12/14 taper)	Quadra-H (12/14 taper)	Quadra-S (12/14 taper)	Quadra-R (12/14 taper)	X-Acta (12/14 taper)
Heads	Ceramtec BIOLOX OPTION (12/14 taper) S Ø 28 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Ceramtec BIOLOX OPTION (12/14 taper) M Ø 28 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Ceramtec BIOLOX OPTION (12/14 taper) L Ø 28 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗
	Ceramtec BIOLOX OPTION (12/14 taper) XL Ø 28 mm	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
	Mectacer BIOLOX OPTION (ball head Ø 28 mm) + sleeve (12/14 taper) size S	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Mectacer BIOLOX OPTION (ball head Ø 28 mm) + sleeve (12/14 taper) size M	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Mectacer BIOLOX OPTION (ball head Ø 28 mm) + sleeve (12/14 taper) size L	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗
	Mectacer BIOLOX OPTION (ball head Ø 28 mm) + sleeve (12/14 taper) size XL	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗

- ✓ - Indicates «Combination approved by Medacta International»
- ✗ - Indicates «Combination NOT approved by Medacta International»
- iiii* only for Mpac Rim-hole acetabular shell from Ø56 to Ø70mm

		Liners							
		Versafitcup Double Mobility UHMWPE/DM Liner UHMWPE	Versafitcup Double Mobility HC UHMWPE/DM Liner HC UHMWPE	Versafitcup CC Flat and Hooded UHMWPE	Versafitcup CC Flat and Hooded HC UHMWPE	Mectacer BIOLOX delta	Ceramtec BIOLOX delta	Mpac Flat and Hooded HC UHMWPE	Double Mobility Converter Liner
Acetabular Shells	Mpac No-hole	✗	✗	✗	✗	✓	✓	✓	✓
	Mpac Two-hole	✗	✗	✗	✗	✓	✓	✓	✓
	Mpac DM	✓	✓	✗	✗	✗	✗	✗	✗
	Mpac Multi-hole	✗	✗	✗	✗	✗	✗	✓	✓
	Mpac Rim-hole	✗	✗	✗	✗	✓ <i>iiii</i>	✓	✓	✓
	Versacem	✓	✓	✗	✗	✗	✗	✗	✗
	Versafitcup CC Trio	✗	✗	✓	✓	✓	✓	✗	✗
	Versafitcup CC Trio No-Hole	✗	✗	✓	✓	✓	✓	✗	✗
	Versafitcup DM	✓	✓	✗	✗	✗	✗	✗	✗